


PRISME TONOMÈTRE À USAGE UNIQUE SINGLE USE TONOMETER PRISM

"Tonoclean" Avec "Clip Insertion"



Sécurité :

Élimine le risque d'infection croisée entre les patients et les protège des maladies transmises par l'intermédiaire du film lacrymal.

Facilité d'utilisation :

Grace au système de clip insertion, le Tonoclean s'installe et se désinstalle aisément sur tous les supports, son conditionnement permet son montage sans contact direct.

Type :

Jetables, stériles, emballés individuellement
Montage : convient au tonomètre type Goldman

Safe

Eliminates the risk of cross-infection between patients and protects them against crossinfection of diseases transmitted through the tear film.

Ease of use

Due to the clip insertion system, the Tonoclean is easy to install and uninstall on all supports, its packaging allows it to be mounted without direct contact.

Type:

Disposable, sterile, individually wrapped
Mounting : suitable for Goldman tonometer type

Conditionnement Packaging	Emballées individuellement sous double sachet stérile par boîte de 5 Individually wrapped in a double sterile bag in box of 5
Stérilisation Sterilization	Ethylene Oxyde

MVCT-01 v2023-06